

An advanced FCR system offering high resolution of 50 μm . Compact to fit anywhere.



Main Features

- 1 Ultra High-resolution image quality
 - 50 micron pitch IP reader
 - Newly developed IP "UR-1"
 - Fuji unique image processing technology
- 2 Cycle time as fast as 97sec. when using 35.4x43cm IP "UR-1".
- 3 15 x 30cm IP "ST-VI" can be read. (*)
- 4 Compact body size saves working space.

(*) reading pitch is 100 micron.

Specifications

Compatible IP	UR-1	35.4x43cm, 24x30cm, 18x24cm
	ST-VI	35.4x43cm/14x7 inch, 24x30cm/10x12inch, 18x24cm/8x10inch, 15x30cm
Reading resolution	50 micron meter (when using UR-1) 100 micron meter (when using ST-VI)	
Cycle time for IP feeding/reading	Approx. 97sec. (35.4x43cm IP UR-1) Approx. 58sec. (35.4x43cm/14x7inch IP ST-VI)	
Reading Grey Scale	12bits/pixel	
Dimensions (W x D x H)	590 x 380 x 810 mm	
Weight	approx. 100kg	

Specifications are subject to change without notice.